# PROVOX<sup>®</sup>

## Measure & Measure Flanges

**Rx  
ONLY**


**UK  
CA**




Measure Flange







Measure

Measure Flange

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner.

The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## CONTENTS

EN - ENGLISH.....	7
DE - DEUTSCH.....	10
NL - NEDERLANDS.....	14
FR - FRANÇAIS.....	17
IT - ITALIANO .....	20
ES - ESPAÑOL.....	24
PT - PORTUGUÊS.....	27
SV - SVENSKA .....	30
DA - DANSK.....	33
NO - NORSK.....	36
FI - SUOMI.....	39
IS - ÍSLENSKA .....	42
LT - LIETUVIŲ K. ....	45
CS - ČESKY .....	48
HU - MAGYAR .....	51
PL - POLSKI.....	54
EL - ΕΛΛΗΝΙΚΑ .....	57
TR - TÜRKÇE.....	61
RU - РУССКИЙ.....	64
JA - 日本語.....	68
KO - 한국어 .....	71
AR - العربية.....	76

## **Intended use**

Provox Measure is intended for sizing the length (corresponding to voice prosthesis length) of tracheoesophageal (TE) punctures.

## **Contraindications**

- Do not use the device on punctures of diameter less than 20 Fr as this may cause damage and/or bleeding of the puncture.
- The device is not intended to be used at the time of surgical creation of the puncture.

## **Device description**

Provox Measure is made of a reusable stainless steel rod and a reusable plastic tube (polyoxymethylene). A single-use silicone flange is attached to the stainless steel rod prior to use.

## **WARNINGS**

- To avoid cross contamination, use appropriate hygienic procedures and make sure that each used device is cleaned, disinfected, dried and steam sterilized before use.
- Disassemble the stainless steel rod and plastic tube prior to cleaning.
- The measure flange is for single use. Re-use and re-processing may cause cross-contamination and damage to the device, which could cause patient harm.
- If the plastic tube does not slide smoothly over the stainless steel rod the complete instrument must be discarded.
- Do not use lubricating (anaesthetizing) gel if the patient has allergies related to these substances.
- Never use a damaged product.

## **PRECAUTIONS**

- Never use excessive force to insert or remove the Measure stainless steel rod and plastic tube as this may cause irritation or bleeding.
- Patients on anticoagulant therapy should be carefully evaluated for the risk of hemorrhage prior to sizing of the TE puncture.
- If lubrication of the flanges is needed, only use water soluble lubricants and make sure that only a thin layer is used.
- Position Provox Measure Flanges correctly according to Figures A and B. Wrong positioning of flanges may result in damage and/or bleeding of the puncture.
- In order to confirm the correct reading, proper fit and leakage around the voice prosthesis shall be checked after insertion of voice prosthesis as recommended in the IFU of the voice prosthesis.

## **Instructions for use**

### **Operating instruction**

Before use always check that the plastic tube slides smoothly over the stainless steel rod. If not, discard the complete instrument.

A single-use silicone flange is attached to the stainless steel rod. The flange shall be positioned according to Figures A and B, respectively depending on if the puncture is 20 Fr or 21-24 Fr. The sizing procedure shall be performed according to Figures 1-6. In order to get a correct reading, push the plastic tube until it gently touches the tracheal mucosa. Note that patients with small stomas might experience high breathing resistance for the duration of the measuring procedure. In order to confirm the correct reading, proper fit and leakage around the voice prosthesis shall be checked after insertion of voice prosthesis as recommended in the IFU of the voice prosthesis.

### **Cleaning and sterilization**

Remove and discard the measure flange after use. Disassemble the stainless steel rod and plastic tube prior

to cleaning. Reassemble the rod and the tube prior to sterilization.

Cleaning, disinfection and sterilization are described in “Instructions for cleaning and sterilization”. If this information is missing, please contact your supplier or visit our website on [www.atosmedical.com](http://www.atosmedical.com).

**CAUTION:** Always ensure mechanical integrity before each use.

### **How supplied**

The device is manufactured in a cleanroom and supplied non-sterile.

### **Device lifetime**

Instruments of stainless steel and polyoxymethylene generally have a long service life. Always ensure mechanical integrity before each use. Provox Measure Flanges are for single use.

### **Disposal**

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

### **Date of printing**

See version number on the back cover of this manual.

### **User assistance information**

For additional help or information please see back cover of this manual for contact information.

### **Reporting**

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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© Coloplast A/S. All rights reserved.
Atos and the Atos Medical logo are trademarks of Coloplast A/S.
**Atos Medical AB**